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Patient Materials

for ophthalmology professionals

- ✓ Educate your patients and save time on consultation
- ✓ Manage patient expectations
- ✓ Successfully grow your business

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DECENTERED IOL - HELP
JUL 17, 2020
A state of ophthalmic readiness: perspectives & actions to return to surgery.

Welcome to the ZEISS webinar:
"Latest clinical data and research on trifocal IOL technology" *
JUN 26, 2020
Latest clinical data and research on trifocal IOL technology

Welcome to the ZEISS webinar:
Mastering Multifocal IOLs "Transition from Cataract Surgery to Cataract Refractive Surgery" *
JUN 26, 2020
Transition from cataract surgery to cataract refractive surgery

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- 1) **Kohnen, T. et al. (2016):** „Trifocal Intraocular Lens Implantation to Treat Visual Demands in Various Distances Following Lens Removal“. *Am J Ophthalmol*
- 2) **Fernández, J. et al. (2021):** „Long-Term Efficacy, Visual Performance and Patient Reported Outcomes with a Trifocal Intraocular Lens: A Six-Year Follow-up“. *Journal of Clinical Medicine*.
- 3) **Boehm, J. et al. (2018):** „Defocus curves of 4 presbyopia-correcting IOL designs: Diffractive panfocal, diffractive trifocal, segmental refractive, and extended-depth-of-focus“. *Journal of Cataract & Refractive Surgery*
- 4) **Plaza-Puche, A. et al. (2016):** „Analysis of defocus curves of different modern multifocal intraocular lenses“. *European Journal of Ophthalmology*
- 5) **Lubinski, W. et al. (2020):** „Comparison of visual outcomes after implantation of At LISA tri 839 MP and Symphony intraocular lenses“. *International Ophthalmology*
- 6) **Meng, J. et al. (2020):** „Decentration and tilt of plate-haptic multifocal intraocular lens in myopic eyes“. *Eye And Vision*
- 7) **Mojzis, P. et al. (2015):** „Visual outcomes of a new toric trifocal diffractive intraocular lens“. *Journal of Cataract & Refractive Surgery*
- 8) **Piovella, M. et al. (2018):** „Patient outcomes following implantation with a trifocal toric IOL: twelve-month prospective multicentre study“. *Eye*

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Publikation	UDVA	UIVA	UNVA	Patienten-zufriedenheit	Unabhängigkeit von Sehhilfen
Alfonso JF, Fernández-Vega Cueto L, Belda-Salmerón L, et al. Visual function after implantation of a diffractive aspheric trifocal intraocular lens. Eur J Ophthalmol 2016;26(5):405–11.	x	x		x	
Alió JL, Kaymak H, Breyer D, et al. Quality of life related variables measured for three multifocal diffractive intraocular lenses: A prospective randomised clinical trial. Clinical & Experimental Ophthalmology 2017.	x	x	x	x	
Al-Khateeb G, Shajari M, Kohnen T. Intraindividual comparative analysis of the visual performance after cataract surgery with implantation of a trifocal and a bifocal intraocular lens. J Cataract Refract Surg 2017;43(5):695–98.	x	x	x	x	
Almulhim AK, Alarfaj KM, Altaisan AA, Alromaih AZ, Aldawod RA. Visual outcomes and patient satisfaction after bilateral implantation of a new trifocal diffractive intraocular lens. Saudi J Ophthalmol. 2018 Oct-Dec;32(4):310-317.	x	x	x	x	x
Ang RE. Comparison of tolerance to induced astigmatism in pseudophakic eyes implanted with small aperture, trifocal, or monofocal intraocular lenses. Clin Ophthalmol. 2019 May 30;13:905-911. doi: 10.2147/OPHTH.S208651. eCollection 2019.	x				
Bilbao-Calabuig R, Llovet-Rausell A, Ortega-Usobiaga J, et al. Visual Outcomes Following Bilateral Implantation of Two Diffractive Trifocal IOLs in 10 084 eyes. American Journal of Ophthalmology 2017.	x	x	x	x	x
Böhm M, Hemkeppler E, Herzog M, Schönbrunn S, de'Lenzo N, Petermann K, Kohnen T. Comparison of a panfocal and trifocal diffractive intraocular lens after femtosecond laser-assisted lens surgery. J Cataract Refract Surg. 2018 Dec;44(12):1454-1462.	x	x	x	x	x
Brito P, Salgado-Borges J, Neves H et al. Light-distortion analysis as a possible indicator of visual quality after refractive lens exchange with diffractive multifocal intraocular lenses. J Cataract Refract Surg 2015;41(3):613–22.	x				

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Publikation	UDVA	UIVA	UNVA	Patienten-zufriedenheit	Unabhängigkeit von Sehhilfen
Chen S-Y, Xie C, Wang Y et al. Full-vision maintenance in extra-high myopia from implantable collamer lens to trifocal intraocular lens implantation. Int J Ophthalmol 2018;11(7):1239–42.		x	x		
Chow SSW, Chan TCY, Ng ALK, Kwok AKH. Outcomes of presbyopia-correcting intraocular lenses after laser in situ keratomileusis. Int Ophthalmol. 2019 May;39(5):1199-1204.	x		x		x
Farideh D, Azad S, Feizollah N et al. Clinical outcomes of new toric trifocal diffractive intraocular lens in patients with cataract and stable keratoconus: Six months follow-up. Medicine (Baltimore) 2017;96(12):e6340.	x	x	x		
Ganesh S, Brar S, Pawar A. Long-term visual outcomes and patient satisfaction following bilateral implantation of trifocal intraocular lenses. Clin Ophthalmol 2017;11:1453–59. 28860693.	x		x	x	x
Gerl M, Breyer DRH, Hagen P et al. Klinischer Vergleich einer trifokalen und einer trifokal-torischen Intraokularlinse auf gleicher diffraktiver Plattform: Clinical Comparison of a Trifocal and a Trifocal-Toric Intraocular Lens Based on the Same Diffractive Platform. Klin Monbl Augenheilkd 2017.	x	x	x		
Höhn F, Tandogan T, Breyer DRH et al. Funktionelle Ergebnisse 1 Jahr nach Implantation einer bitorischen, trifokalen Intraokularlinse. Klin Monbl Augenheilkd 2015;232(8):957–61.	x	x	x		
Kim BH, Hyon JY, Kim MK. Effects of Bifocal versus Trifocal Diffractive Intraocular Lens Implantation on Visual Quality after Cataract Surgery. Korean J Ophthalmol. 2019 Aug;33(4):333-342.	x	x	x	x	
Kim M, Kim JH, Lim TH et al. Comparison of Reading Speed after Bilateral Bifocal and Trifocal Intraocular Lens Implantation. Korean J Ophthalmol 2018;32(2):77–82.	x	x	x		
Kohnen T, Titke C, Böhm M. Trifocal Intraocular Lens Implantation to Treat Visual Demands in Various Distances Following Lens Removal. American Journal of Ophthalmology 2016;161:71-7.e1 (accessed 19 Apr 2016).	x	x	x	x	x

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Publikation	UDVA	UIVA	UNVA	Patienten-zufriedenheit	Unabhängigkeit von Sehhilfen
Kretz FTA, Choi CY, Müller M et al. Visual Outcomes, Patient Satisfaction and Spectacle Independence with a Trifocal Diffractive Intraocular Lens. Korean J Ophthalmol 2016;30(3):180–91.	x	x	x	x	x
Kretz, Florian T A, Breyer D, Diakonis VF et al. Clinical Outcomes after Binocular Implantation of a New Trifocal Diffractive Intraocular Lens. Journal of Ophthalmology 2015;2015:962891	x	x	x		
Kretz, Florian T A, Breyer D, Klabe K et al. Clinical Outcomes After Implantation of a Trifocal Toric Intraocular Lens. J Refract Surg 2015;31(8):504–10	x	x	x	x	x
Kretz, Florian T A, Müller M, Gerl M et al. Binocular function to increase visual outcome in patients implanted with a diffractive trifocal intraocular lens. BMC Ophthalmol 2015;15:110.	x	x	x		
Law EM, Aggarwal RK, Kasaby H. Clinical outcomes with a new trifocal intraocular lens. Eur J Ophthalmol 2014;24(4):501–08	x	x	x	x	
Lenton L. Visual performance in a flight simulator: multifocal intraocular lenses in pilots. BMJ Open Ophthalmol 2018;3(1):e000139.	x	x	x		x
Liu X, Xie L, Huang Y. Comparison of the Visual Performance After Implantation of Bifocal and Trifocal Intraocular Lenses Having an Identical Platform. J Refract Surg 2018;34(4):273–80.		x	x	x	
Marques EF, Ferreira TB. Comparison of visual outcomes of 2 diffractive trifocal intraocular lenses. J Cataract Refract Surg 2015;41(2):354–63. 25661129	x	x	x		x
Martínez de Carneros-Llorente A, Martínez de Carneros A, Martínez de Carneros-Llorente P, Jiménez-Alfaro I. Comparison of visual quality and subjective outcomes among 3 trifocal intraocular lenses and 1 bifocal intraocular lens. J Cataract Refract Surg. 2019 May;45(5):587-594.	x			x	
Martínez-de-la-Casa JM, Carballo-Alvarez J, García-Bella J et al. Photopic and mesopic performance of 2 different trifocal diffractive intraocular lenses. Eur J Ophthalmol 2016:0.	x	x	x		

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Publikation	UDVA	UIVA	UNVA	Patienten-zufriedenheit	Unabhängigkeit von Sehhilfen
Mencucci R, Favuzza E, Caporossi O et al. Comparative analysis of visual outcomes, reading skills, contrast sensitivity, and patient satisfaction with two models of trifocal diffractive intraocular lenses and an extended range of vision intraocular lens. Graefes Arch Clin Exp Ophthalmol 2018.	x	x	x	x	x
Mencucci R, Favuzza E, Caporossi O et al. Visual performance, reading ability and patient satisfaction after implantation of a diffractive trifocal intraocular lens. Clin Ophthalmol 2017;11:1987–93.	x	x	x	x	x
Mendicute J, Kapp A, Levy P et al. Evaluation of visual outcomes and patient satisfaction after implantation of a diffractive trifocal intraocular lens. J Cataract Refract Surg 2016;42(2):203–10.				x	x
Mojzis P, Kukuckova L, Majerova K et al. (2017). Postoperative visual performance with a bifocal and trifocal diffractive intraocular lens during a 1-year follow-up. Int J Ophthalmol, 17 October 2017:1528–33.	x	x	x		
Mojzis P, Kukuckova L, Majerova K et al. Comparative analysis of the visual performance after cataract surgery with implantation of a bifocal or trifocal diffractive IOL. J Refract Surg 2014;30(10):666–72.	x	x	x	x	
Mojzis P, Majerova K, Hrcckova L et al. Implantation of a diffractive trifocal intraocular lens: one-year follow-up. J Cataract Refract Surg 2015;41(8):1623–30	x	x	x	x	
Mojzis P, Majerova K, Plaza-Puche AB et al. Visual outcomes of a new toric trifocal diffractive intraocular lens. J Cataract Refract Surg 2015;41(12):2695–706.	x		x		
Mojzis P, Peña-García P, Liehneova I et al. Outcomes of a new diffractive trifocal intraocular lens. J Cataract Refract Surg 2014;40(1):60–69	x	x	x		
Piovela M, Colonval S, Kapp A et al. Patient outcomes following implantation with a trifocal toric IOL: twelve-month prospective multicentre study. Eye (Lond)				x	x
Plaza-Puche AB, Alio JL, Sala E et al. Impact of low mesopic contrast sensitivity outcomes in different types of modern multifocal intraocular lenses. Eur J Ophthalmol 2016:0.	x		x		

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Publikation	UDVA	UIVA	UNVA	Patienten-zufriedenheit	Unabhängigkeit von Sehhilfen
Plaza-Puche AB, Alio JL. Analysis of defocus curves of different modern multifocal intraocular lenses. Eur J Ophthalmol 2016;0.	x				
Postolache C, Postolache O. COMPARATION OF REFRACTIVE RESULTS WITH BIFOCAL IMPLANTS AT LISA 809 AND TRIFOCAL AT LISA TRI839. Rom J Ophthalmol 2015;59(2):100–02	x	x	x		
Qi Y, Lin J, Leng L et al. Role of angle κ in visual quality in patients with a trifocal diffractive intraocular lens. J Cataract Refract Surg 2018;44(8):949–54.		x	x	x	
Sezgin Asena B. Visual and refractive outcomes, spectacle independence, and visual disturbances after cataract or refractive lens exchange surgery: Comparison of 2 trifocal intraocular lenses. J Cataract Refract Surg. 2019 Nov;45(11):1539-1546.	x	x	x		
Steinwender G, Schwarz L, Böhm M et al. Visual results after implantation of a trifocal intraocular lens in high myopes. J Cataract Refract Surg 2018;44(6):680–85.	x	x	x		
Tarib I, Diakonov VF, Breyer D, Höhn F, Hahn U, Kretz FTA. Outcomes of combining a trifocal and a low-addition bifocal intraocular lens in patients seeking spectacle independence at all distances. J Cataract Refract Surg. 2019 May;45(5):620-629.	x	x	x		x
Velasco-Barona C, Corredor-Ortega C, Mendez-Leon A, Casillas-Chavarin N L, Valdepeña-Lopez Velarde D, Cervantes-Coste G, Malacara-Hernandez D und Gonzalez-Salinas R. Influence of angle κ and higher-order aberrations on visual quality employing two diffractive trifocal IOLs: (AT LISA tri 839MP vs Acrysof IQ PanOptix). Journal of Ophthalmology 2019, Article ID 7018937	x	x	x		
Wang W, Ni S, Li X et al. Femtosecond laser-assisted cataract surgery with implantation of a diffractive trifocal intraocular lens after laser in situ keratomileusis: a case report. BMC Ophthalmol 2018;18(1):160.	x	x	x	x	
Xue S, Zhao G, Yin X et al. Effect of incision on visual outcomes after implantation of a trifocal diffractive IOL. BMC Ophthalmol 2018;18(1):171.	x	x	x		
Yang Y, Lv H, Wang Y et al. Clinical outcomes following trifocal diffractive intraocular lens implantation for age-related cataract in China. Clin Ophthalmol 2018;12:1317–24.	x	x	x	x	x

Seeing beyond