


FORM FIT[®]
INTRACANALICULAR PLUG

A New Advance in Dry Eye Treatment




ONE SIZE. ONE SOLUTION.

OASIS[®] 

The future in punctal occlusion is here.

OASIS introduces the next generation of punctal occlusion... Form Fit®.



INSERTION OF THE FORM FIT®

Form Fit® is easily inserted into the canaliculus. Form Fit® comes pre-loaded on a user-friendly inserter. The polyimide sheath is inserted into the punctum and the inserter handle compressed to expel the hydrogel plug. Once the hydrogel plug is expelled, maintain pressure on the inserter handle to remove the polyimide sheath from the punctum.

In its dry state, the plug measures 0.3mm in diameter by 3.0mm long. Form Fit® easily slides into the punctal cavity. The need for punctal measurement has been eliminated. Form Fit® then swells to form a gelatinous material.



ONE SIZE. ONE SOLUTION.

Form Fit® is truly a one size plug solution. Once inserted into the canaliculus, Form Fit® hydrates over a 10 minute period. As it hydrates, Form Fit® increases in size until it completely fills the vertical canalicular cavity to form a perfect fit.

HOW IT WORKS

Form Fit® is made of a hydrogel material that expands into a soft, pliable, gelatinous material when it contacts tear film. It takes approximately 10 minutes for the plug to become completely hydrated.

Form Fit® securely conforms to the punctal cavity so there is no need for a cap or anchor on the end of the plug. This greatly increases patient comfort and decreases foreign body sensation.

Form Fit® does not begin to swell until it makes contact with tear film reducing the potential for product waste. Form Fit® is made from a hydrogel material that is safe and biocompatible.

REMOVAL OF THE FORM FIT®

To remove Form Fit®, simply flush saline solution through the punctal opening and the gelatinous material will be expressed through either the nasal cavity or down the throat.

BUSINESS SENSE

- Form Fit® requires no special storage and reacts consistently every time.
- Plug inventory is reduced as multiple plug sizes are no longer necessary.
- Removing the need for measurement provides quicker exam time.
- Form Fit® expands beyond other intracanalicular plug diameters

PRODUCT FEATURES

- Hydrates over a period of 10 minutes
- Expands to softly fit the individual's punctal cavity without overstretching the tissue
- Eliminates the cap and anchor system traditionally seen (silicone plugs)
- Easily removed via irrigation of the tear duct
- Biocompatible
- No foreign body sensation
- Conveniently pre-loaded on a user-friendly inserter
- Sizing of punctum is eliminated



OASIS offers a complete solution to punctal occlusion with the **SOFT PLUG®** family of **punctum plugs**.



COLLAGEN INTRACANALICULAR PLUGS

The collagen plugs have a smooth proprietary finish that simplifies insertion of the plug. Collagen plugs are ideal for short-term, diagnostic and post-surgical occlusion. Diameter Sizes: 0.2mm, 0.3mm and 0.4mm



EXTENDED DURATION PLUGS

Extended Duration absorbable plugs are designed to provide patient comfort for up to 3 months. This long lasting plug is effective and economical for transient dry eye symptoms following surgery. Diameter Sizes: 0.3mm, 0.4mm and 0.5mm



SILICONE PUNCTUM PLUGS

SOFT PLUG® Silicone plugs feature a softer more flexible silicone that maximizes patient comfort. Sizes: Micro (0.4mm), Mini (0.5mm), Petite (0.6mm), Small (0.7mm) and Medium (0.8mm)

To Order Call 909.305.5400 or 800.528.9786

24-hour FAX 909.305.9987 or 800.631.7210

www.OasisMedical.com

514 S. Vermont Avenue
Glendora, CA 91741 U.S.A.

