



Preloaded Hydrophobic Soft Acrylic



Lens Model	N4-11B
Material Optic	Hydrophobic soft acrylic resin
Material Haptic	Polymethylmethacrylate (PMMA)
Optic Design	Biconvex
Optic Length	6.0 mm
Overall Diameter	12.5 mm
Haptic Angle	7°
A-Constant*	119.2
Diopter Range	+10.0~27.0D (0.5D increments) +27.0~28.0D (1.0D increments)

IOL Constant for Zeiss IOL-Master					
*Nominal	Haigis	Hoffer Q	Holl. 1	SRK/T	SRKII
A=119.2	a0=1.8386 a1=0.40 a2=0.10	pACD=6.0312	sf=2.2492	A=119.8	A=120.4

Manufactured by:
STAAR Japan Inc., 717-30, Futamata, Ichikawa-shi, Chiba 272-0001, Japan

EU Authorized Representative:
STAAR Surgical AG, Hauptstrasse 104, PO Box 463, CH-2560 Nidau, Switzerland
Phone +41 32 332 88 88. Fax +41 32 332 88 99, info@staarag.ch www.staar.com

