




Datenblatt Kapselspannringe

Lenstec offers Capsular Tension Rings in both clear and blue PMMA. They offer high flexibility, high memory control and are easy to insert manually or with any available injector.

Both designs are applicable when the zonules are missing or damaged; in cases of lens luxation during Phaco complications; zonular breakage during sudden unfolding of a foldable lens; for repositioning of a lens or if decentration of same. Other uses apply for myopia, pseudoexfoliation, eyes undergoing vitrectomy and trauma.

The implantation of each design offers various advantages: force distribution acting circumferentially on the capsule; less risk of capsular shrinking; less risk of capsular fibrosis; and safer IOL centration



Model No.	Size		Color
	Open	Compressed	
LR-1300	12.30	10.00	Clear
LR-1400	14.30	12.00	Clear
LR-1300B	12.30	10.00	Blue
LR-1400B	14.30	12.00	Blue

LR-1400 and LR-1400B are used in highly myopic eyes (IOL power <+22D)