

# Axial Marker

DESIGNED BY P. FAZIO M.D.



## J 2281.2

Manual axial marker for horizontal axis (0°/180°),  
for patient seated at the slit-lamp.  
(based on Elies M.D. model)

Allows two preoperative marking points (0°/180°)

Curve shaped for better visibility during the marking procedure

Certified  
Quality System

ISO 9001  
ISO 13485

© Copyright:  
Hermann Janach



E. Janach S.r.l.  
Via Borgo Vico, 35 - 22100 Como - Italy  
Tel. +39 031 574088 - fax +39 031 572055  
international@janach.it - [www.janach.it](http://www.janach.it)

**e.janach®**

