
















# OASIS Surgical

## NEW DISPOSABLE INSTRUMENTS

- Packaged sterile 10 per box or bulk
- Also available non-sterile in bulk

		REF #	GAUGE	LENGTH	DESCRIPTION
ANESTHESIA		4004H	27g	1.25"/32mm	Atkinson Tip
		4004K	30g	0.625"/16mm	Sharp Tip
VISCO CANNULAS		4003D	20g	0.875"/22mm	Curved
		4002H	25g	0.875"/22mm	Angled 8mm from tip
		4002D	27g	0.875"/22mm	Angled 10mm from tip
IRRIGATION		4001D	21g	0.875"/22mm	Angled 4mm from tip
		4002E	27g	0.875"/22mm	Angled 5mm from tip
		4002R	27g	0.875"/22mm	Formed 11mm from end
		4002S	27g	0.875"/22mm	Straight
CAPSULOTOMY		4017K	27g	0.5"/13mm	10.5mm bend from tip
HYDRO- DISSECTION		4068D	27g	0.875"/22mm	Curved
CAPSULE POLISHERS		4027F	23g	0.875"/22mm	Curved; 0.33mm end port Tip: 2mm long x 1mm wide
		4028D	27g	0.875"/22mm	Curved; 23mm soft tip, beveled down
LACRIMAL		6010	26g	1.125"/28mm	Curved; open tip
VITREORETINAL		4073	20g	1.125"/28mm	Straight

DISTRIBUTED BY:

Manufactured by:  
OASIS Medical, Inc.  
Glendora, California  
Tel: 909.305-5400 ■ Fax: 909.305.9987  
www.oasismedical.com