

OASIS[®] Surgical

DISPOSABLE SURGICAL PRODUCTS



Introducing the Premier Edge® Safety Knives



Patent # D537528

Designed for safety and efficiency

- ▶ Transparent safety guard slides fully over blade when not in use
- ▶ Fully retracts when blade is in use
- ▶ Extend, retract, and lock the guard in place with one hand
- ▶ Clear View of the surgical field
- ▶ Non-slip grip
- ▶ Sterile; 6 knives per box
- ▶ Bulk and non-sterile packaging also available



REF #	DESCRIPTION	SIZE
PE3015-SL	Safety Stab Knife	15°
PE3520-SL	Safety MVR Knife	20g
PE3720A-SL	Safety Crescent Knife, Angled	2.0mm
PE4822-SL	Safety Thin Line Slit Knife	2.2mm
PE4824-SL	Safety Thin Line Slit Knife	2.4mm
PE4826-SL	Safety Thin Line Slit Knife	2.65mm
PE4827-SL	Safety Thin Line Slit Knife	2.75mm
PE4828-SL	Safety Thin Line Slit Knife	2.85mm
PE4830-SL	Safety Thin Line Slit Knife	3.0mm
PE4832-SL	Safety Thin Line Slit Knife	3.2mm

OASIS
Surgical

www.oasismedical.com

Premier Edge® Microsurgical Knives - Disposable

Ergonomic Handle Designed for comfort, safety, efficiency and ease of use to provide the most natural feel in the surgeon's hand and enhance overall control during surgery

Blade Consistently sharp with non-glare finish creates reproducible results

Packaging Sterile blister packs, 6 per box

STAB



**SLIT
THIN LINE**



CLEAR CORNEA



PARALLEL SIDE



TRAPEZOID



IMPLANT



SIDE PORT



ROUND TUNNEL



CRESCENT TUNNEL



MVR






REF #	DESCRIPTION	WIDTH	ANGLE	BEVEL
STAB KNIVES				
PE3015	Stab	15°		Double
PE3015-3	Stab	15°	3mm depth	
PE3015-5	Stab	15°	5mm depth	
PE3022	Stab	22.5°		Double
PE3030	Stab	30°		Double
PE3045	Stab	45°		Double
THIN LINE SLIT KNIVES				
PE4822-TL	Thin Line Slit	2.2mm	Angled	Bevel Up
PE4824-TL	Thin Line Slit	2.4mm	Angled	Bevel Up
PE4826-TL	Thin Line Slit	2.65mm	Angled	Bevel Up
PE4827-TL	Thin Line Slit	2.75mm	Angled	Bevel Up
PE4828-TL	Thin Line Slit	2.85mm	Angled	Bevel Up
PE4830-TL	Thin Line Slit	3.0mm	Angled	Bevel Up
PE4832-TL	Thin Line Slit	3.2mm	Angled	Bevel Up
CLEAR CORNEA KNIVES				
PE3122	Clear Cornea	2.2mm	Angled	Double
PE3125	Clear Cornea	2.5mm	Angled	Double
PE3126	Clear Cornea	2.65mm	Angled	Double
PE3127	Clear Cornea	2.75mm	Angled	Double
PE3128	Clear Cornea	2.85mm	Angled	Double
PE3132	Clear Cornea	3.2mm	Angled	Double
CLEAR CORNEA PARALLEL SIDE				
PE3226	Parallel Side	2.65mm	Angled	Double
PE3227	Parallel Side	2.75mm	Angled	Double
PE3228	Parallel Side	2.85mm	Angled	Double
PE3230	Parallel Side	3.0mm	Angled	Double
CLEAR CORNEA TRAPEZOID				
PE3301	Trapezoid	1.8mm – 2.2mm	Angled	Bevel Up
PE3302	Trapezoid	2.5mm – 3.0mm	Angled	Bevel Up
IMPLANT KNIVES				
PE3440	Implant	4.0mm	Angled	Bevel Up
PE3452	Implant	5.2mm	Angled	Bevel Up
SIDE PORT KNIVES				
PE3010	Side Port	1.0mm	Angled	Double
PE3012	Side Port	1.2mm	Angled	Double
PE3014	Side Port	1.4mm	Angled	Double
ROUND TUNNEL KNIVES				
PE3622	Round Tunnel	2.2mm	Angled	Double
PE3630*	Round Tunnel	3.0mm	Angled	Double
CRESCENT TUNNEL KNIVES				
PE3720	Crescent Tunnel	2.0mm	Straight	Bevel Up
PE3720A	Crescent Tunnel	2.0mm	Angled	Bevel Up
MVR KNIVES				
PE3519	MVR	19g	Straight	Double
PE3519A	MVR	19g	Angled	Double
PE3520	MVR	20g	Straight	Double
PE3520A	MVR	20g	Angled	Double
PE3523	MVR	23g	Straight	Double
PE3523A	MVR	23g	Angled	Double
PE3525	MVR	25g	Straight	Double

* Make to Order



Needles, Cannulas, Cystotomes – Disposable




■ Available sterile 10 per box**

Gauge	Anesthesia 	Visco/Irrigation Cannulas 	Capsulotomy Instruments 
19g (1.10)	4005 (1"/25mm) 4005TP (1"/25mm)	4036 (1"/25mm)	
20g (0.90)		* 4003D (0.875"/22mm) Curved	
21g (0.80)		* 4001D (0.875"/22mm) 4mm Bend	
23g (0.60)	4004C (1.5"/38mm) Atkinson 4004F (1.25"/32mm) Atkinson	4001B (0.875"/22mm) 10mm Bend * 4001CF (0.875"/22mm) 4mm Bend 4001CTW (0.875"/22mm) 8mm Thin Wall	
24g (0.55)			
25g (0.50)	4004 (1.5"/38mm) 4004B (1.5"/38mm) Atkinson 4004D (0.875"/22mm) Atkinson 4004M (1.25"/32mm) Atkinson	4001A (0.875"/22mm) 6mm Bend 4001AF (0.875"/22mm) 4mm Bend * 4002H (0.875"/22mm) 8mm Bend 4003 (1.125"/29mm) Curved	4017 (0.875"/22mm) Formed 4017A (0.875"/22mm) Straight 4017G (0.625"/16mm) Formed 4017IE (0.625"/16mm) Formed 4017L (0.625"/16mm) Reverse Formed 4017M (0.625"/16mm) Angled 4017S (0.625"/16mm) Formed 4017Z (0.625"/16mm) F-SL Curve
26g (0.45)			
27g (0.40)	* 4004H (1.25"/32mm) Atkinson	4002 (0.875"/22mm) 6mm Bend 4002A (0.875"/22mm) 11mm Bend 4002B (0.875"/22mm) 8mm Bend * 4002D (0.875"/22mm) 10mm Bend * 4002E (0.875"/22mm) 5mm Bend 4002F (0.875"/22mm) 4mm Bend 4002L (0.875"/22mm) 9mm Bend * 4002R (0.875"/22mm) Formed 11mm * 4002S (0.875"/22mm) Straight	4017B (0.875"/22mm) Formed 4017IS (0.625"/16mm) Straight 4017IO (0.625"/16mm) Formed * 4017K (0.5"/13mm) 10.5mm Bend 4017N (0.625"/16mm) Angled 4017R (0.625"/16mm) Formed * 4017Y (0.625"/16mm) F-SL Curve
30g (0.30)	* 4004K (0.625"/16mm) Sharp Tip	4001 (0.875"/22mm) 6mm Bend 4001F (0.875"/22mm) 4mm Bend	4017I (0.625"/16mm) Formed 4017IR (0.625"/16mm) Reverse Formed

* Make to Order

Cannulas, Cortex Extractors, Capsule Polishers – Disposable

■ Available sterile 10 per box*

Gauge	Hydrodissection/LASIK Cannulas	Cortex Extractors and Vectis	Capsule Polishers
19g (1.10)			
20g (0.90)			
21g (0.80)		4016 (0.875"/22mm) U-shaped	4027 (0.875"/22mm) 9mm Bend 4027B (0.875"/22mm) 6mm Bend
23g (0.60)	4123 (0.875"/22mm)	4007 Simcoe Double I and A 4008 (0.5"/13mm) Simcoe 4016F (0.875"/22mm) U-shaped 4026 (0.5"/13mm) Charleux	4027A (0.875"/22mm) 9mm Bend * 4027F (0.875"/22mm) Olive Tip 4028B (0.875"/22mm) Silicone Sleeve
24g (0.55)		4029A Vectis, three 0.25mm ports	
25g (0.50)	4016AH (0.875"/22mm) U-shaped 4036A (0.875"/22mm) 8mm Bend 4036C (0.875"/22mm) 4mm Bend 4036F (0.875"/22mm) Helsinki 4036M (0.875"/22mm) Sauter 4068 (0.875"/22mm) Kellan 4036L (0.875"/22mm) Lin LASIK 4036LB (0.875"/22mm) Buratto LASIK 4036LR (0.875"/22mm) Lindstrom LASIK 4225 (0.875"/22mm) LASIK, 6 ports	4016G (0.875"/22mm) J-shaped	4027E (1"/25mm) Olive Tip
26g (0.45)	4036LS (1"/25mm) Slade LASIK		
27g (0.40)	4036B (0.875"/22mm) 8mm Bend 4036D (0.875"/22mm) 4mm Bend 4036G (0.875"/22mm) Helsinki 4036N (0.875"/22mm) Sauter 4036J (0.75"/19mm) Chang 4036P (0.875"/22mm) Cleaving, Round 4036Q (0.875"/22mm) Cleaving, Flat * 4068D (0.875"/22mm) Curved	4016B (0.875"/22mm) U-shaped	4028 (0.875"/22mm) Terry * 4028D (0.875"/22mm) Curved, Soft Tip
30g (0.30)			

** Also available in the U.S. through the Bulk Program


Visco Shield® HPMC Viscoelastic

- Available in 20,000 and 8,000 cps
- Stores at room temperature
- Sterile, 6 per box with 23g x 7/8 thin wall cannula

Visco Shield® TOPICAL

- Designed for topical use to hydrate and protect the outer surface of the eye
- Sterile, 6 per box with 23g x 7/8 thin wall cannula


REF#	DESCRIPTION
Viscoelastic – 2ml syringe	
5122	Visco Shield Viscoelastic 20,000 cps
51082	Visco Shield Viscoelastic 8,000 cps
Topical – 1.5ml syringe	
5200	Visco Shield Topical 8,000 cps







Iris Retractors – Disposable

- To temporarily stretch and retain the iris during surgery
- Flexible blue polypropylene hooks
- Flat silicone stopper adjusts easily to secure hook

REF #	DESCRIPTION
4465	Iris Retractors Sterile, 30 per box



Accessories – Disposable

	REF #	DESCRIPTION	QUANTITY
Corneal Marking Pen			
	0442	Corneal Marking Pen, Fine Tip, Gentian Violet Sterile.	10 per box
Disc Fluid Filter			
	0601	For use between syringe and needle. Filters out particulate and biological contaminate matter larger than 0.22 µ. Bidirectional filter: female luer lock and male luer slip with screw-in lock to secure cannula in place. Sterile.	50 per box
Eye Cautery			
	3501-A	Fine tip, battery operated. Low temperature. Sterile. Range: 1100°F ± 50°F	10 per box
SOFT CELL® PVA Foam			
	0525	Surgical Spears Sterile. High density, enhanced fluid capacity, fiber-free	150 per box
	0522-5	LASIK Eye Drain (Chayet) Sterile, 10mm internal diameter	20 per box
	0523-5	LASIK Eye Drain (Chayet/Katzen) Sterile, 11.5mm internal diameter	20 per box
	0529	Wick + Wipe	20 per box

OASIS® 

Glendora, CA

909.305.5400 ■ 909.305.9987 FAX

international@oasismedical.com

www.oasismedical.com

European Authorized Representative
Donawa Lifescience Consulting Srl
Piazza Albania 10
00153 Rome, Italy



ISO 13485:2003

CM004 10/11