


Introducing

VISCO SHIELD® TOPICAL

Designed for **topical use** to hydrate and protect the outer surface of the eye

Visco Topical Layer



Visco Shield Topical is preservative-free and designed to keep the cornea hydrated for up to 20 minutes

Visco Shield Topical adheres to the cornea and forms a thin, optically neutral clear film



- 8,000 cps **HPMC** viscoelastic
 - Optically clear
- Stores at room temperature
- Packaged sterile in syringe, 1.5ml fill
- complete with cannula for easy application
- 6 per box

OASIS® 

Glendora, CA 91741
800.528.9786 ■ 909.305.5400
FAX 800.631.7210 ■ 909.305-9987
www.oasismedical.com